

Anti-Phospho-NEDD4L (S448) Rabbit Polyclonal Antibody



Catalog #: U1931

Aliases

E3 ubiquitin-protein ligase NEDD4-like (EC 6.3.2.-) (NEDD4.2) (Nedd4-2)

Background

Gene Name: NEDD4L/KIAA0439/NEDL3

NCBI Gene Entry: [23327](#)

UniProt Entry: [Q96PU5](#)

Application Information

Molecular Weight: Predicted, 112 kDa; observed, 110-120 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-NEDD4L (Ser448)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of 1 HeLa treated with LPS, 2 HeLa, using primary antibody at 1:1,000 dilution. Secondary antibody was diluted at 1:10,000