

Anti-Phospho-ROS1 (Y2274) Rabbit Polyclonal Antibody



Catalog #: U1971

Aliases

Proto-oncogene tyrosine-protein kinase ROS (EC 2.7.10.1) (Proto-oncogene c-Ros) (Proto-oncogene c-Ros-1) (Receptor tyrosine kinase c-ros oncogene 1) (c-Ros receptor tyrosine kinase)

Background

Gene Name: ROS1/MCF3/ROS

NCBI Gene Entry: [6098](#)

UniProt Entry: [P08922](#)

Application Information

Molecular Weight: Predicted, 264 kDa; observed, 258 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-ROS1 (Tyr2274)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Immunocytochemistry (IC): 1:50-1:200

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HeLa using primary antibody at 1:1,000 dilution 4°C, overnight.
Secondary antibody was diluted at 1:10,000 25°C □ 1.5hours