

Anti-Phospho-FGFR3 (Y760) Rabbit Polyclonal Antibody



Catalog #: U2132

Aliases

Fibroblast growth factor receptor 3 (FGFR-3) (EC 2.7.10.1) (CD antigen CD333)

Background

Gene Name: FGFR3/JTK4

NCBI Gene Entry: [2261](#)

UniProt Entry: [P22607](#)

Application Information

Molecular Weight: Predicted, 88 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FGFR3 (Tyr760)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SY-SY5Y using primary antibody at 1:1,000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10,000 25°C □ 1.5 hours