

Anti-Phospho-ROCK2 (Y722) Rabbit Polyclonal Antibody



Catalog #: U2164

Aliases

Rho-associated protein kinase 2 (EC 2.7.11.1) (Rho kinase 2) (Rho-associated, coiled-coil-containing protein kinase 2) (Rho-associated, coiled-coil-containing protein kinase II) (ROCK-II) (p164 ROCK-2)

Background

Gene Name: ROCK2/KIAA0619

NCBI Gene Entry: [9475](#)

UniProt Entry: [O75116](#)

Application Information

Molecular Weight: Predicted, 161 kDa; observed, 161 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-ROCK2 (Tyr722)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of mouse spleen /Hela ,using primary antibody at 1:1,000 dilution.
Secondary antibody was diluted at 1:10,000