

Anti-Phospho-PDE3A (S312) Rabbit Polyclonal Antibody



Catalog #: U2187

Aliases

cGMP-inhibited 3',5'-cyclic phosphodiesterase A (EC 3.1.4.17) (Cyclic GMP-inhibited phosphodiesterase A) (CGI-PDE A)

Background

Gene Name: PDE3A

NCBI Gene Entry: [5139](#)

UniProt Entry: [Q14432](#)

Application Information

Molecular Weight: Predicted, 126 kDa; observed, 126 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PDE3A (Ser312)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HeLa using primary antibody at 1:1,000 dilution 4°C, overnight.
Secondary antibody was diluted at 1:10,000 25°C □ 1.5hours