

Anti-Phospho-BCLAF1 (S531) Rabbit Polyclonal Antibody



Catalog #: U2231

Aliases

Bcl-2-associated transcription factor 1 (Btf)

Background

Gene Name: BCLAF1/BTF/KIAA0164

NCBI Gene Entry: [9774](#)

UniProt Entry: [Q9NYF8](#)

Application Information

Molecular Weight: Predicted, 106 kDa; observed, 106 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-BCLAF1 (Ser531)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of various,using primary antibody at 1:1,000 dilution. Secondary antibody was diluted at 1:10,000