

Anti-Alpha-tubulin (Acetyl-K112) Rabbit Polyclonal Antibody



Catalog #: U2258

Aliases

TUBA1A; Tubulin Alpha 1a; TUBA3; B-ALPHA-1; Tubulin, Alpha, Brain-Specific; Tubulin Alpha-1A Chain; Tubulin Alpha-3 Chain; Tubulin B-Alpha-1; FLJ25113; Alpha-Tubulin; Hum-A-Tub1; Hum-A-Tub2; EC 3.6.5.-; LIS3; TUBA1B; Tubulin Alpha 1b; K-ALPHA-1; Tubulin Alpha-Ubiquitous Chain; Tubulin, Alpha, Ubiquitous; Alpha-Tubulin Ubiquitous; Tubulin Alpha-1B Chain; Tubulin K-Alpha-1; Alpha Tubulin; EC 3.6.5.-; TUBA1C; Tubulin Alpha 1c; Bcm948; TUBA6; Tubulin Alpha-1C Chain; Tubulin Alpha-6 Chain; Tubulin, Alpha 6; Alpha-Tubulin 6; MGC14580; MGC10851; EC 3.6.5.-

Background

Gene Name: TUBA1A/TUBA1B/TUBA1C

NCBI Gene Entry: [7846/10376/84790](#)

UniProt Entry: [Q71U36/P68363/Q9BQE3](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 55 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Alpha-tubulin (Acetyl-K112)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-Alpha-tubulin (Acetyl-K112) Rabbit Polyclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-Alpha-tubulin (Acetyl-K112) antibody (Cat#U2258). Total lysates (30 µg) were loaded and separated by SDS-PAGE. The blot was incubated with anti-Alpha-tubulin (Acetyl-K112) antibody (Cat#U2258, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Alpha-tubulin (Acetyl-K112) Rabbit Polyclonal Antibody

PAGE 3

Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM