

Anti-Phospho-Casein Kinase II α (Y255) Rabbit Polyclonal Antibody



Catalog #: U2263

Aliases

CSNK2A1; Casein Kinase 2 Alpha 1; Casein Kinase II Subunit Alpha; Casein Kinase 2, Alpha 1 Polypeptide; EC 2.7.11.1; CK2A1; Cka1; Cka2; Casein Kinase II Alpha 1 Polypeptide Pseudogene; Casein Kinase II Alpha 1 Subunit; CK2 Catalytic Subunit Alpha; Protein Kinase CK2; CK II Alpha; EC 2.7.11; OCNDS; CKII; CKA1; CKA2

Background

Gene Name: CSNK2A1

NCBI Gene Entry: [1457](#)

UniProt Entry: [P68400](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-Casein Kinase II α (Tyr255)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:100-1:300

Immunocytochemistry (IC): 1:50-1:200

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Casein Kinase II α (Y255) Rabbit Polyclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM