

Anti-Phospho-DNA-PKCS (S2056) Rabbit Polyclonal Antibody



Catalog #: U2273

Aliases

PRKDC; Protein Kinase, DNA-Activated, Catalytic Subunit; DNA-PKcs; DNPK1; P460; DNA-Dependent Protein Kinase Catalytic Subunit; DNA-PKC; DNAPKc; DNAPK; XRCC7; HYRC1; P350; HYRC; Protein Kinase, DNA-Activated, Catalytic Polypeptide; DNA-PK Catalytic Subunit; Ser-473 Kinase; EC 2.7.11.1; S473K; Hyper-Radiosensitivity Of Murine Scid Mutation, Complementing 1; IMD26

Background

Gene Name: PRKDC

NCBI Gene Entry: [5591](#)

UniProt Entry: [P78527](#)

Application Information

Molecular Weight: Predicted, 469 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-DNA-PKCS (S2056)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:300

Immunocytochemistry (IC): 1:50-1:200

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-DNA-PKCS (S2056) Rabbit Polyclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM