

Anti-CD69 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: UY0002

Aliases

Very Early Activation Antigen (VEA), Activation inducer molecule (AIM)

Background

Gene Name: CD69

NCBI Gene Entry: [969](#)

UniProt Entry: [Q07108](#)

Application Information

Molecular Weight: Predicted, 23 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7255

Species Reactivity: Human, monkey

Applications Tested: Flow cytometry (FCM)

Conjugation: PE

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2%BSA, and 0.1%Proclin 300

Storage

Store at 4 °C in dark.

Recommended Dilutions

5µl/□

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD69 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM