

Anti-CD45RA Mouse Monoclonal Antibody



Catalog #: UY0003

Aliases

GP180, L-CA, LCA, LY5, T200, PTPRC

Background

Gene Name: PTPRC

NCBI Gene Entry: [5788](#)

UniProt Entry: [P08575](#)

Application Information

Molecular Weight: Predicted, 147 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7260

Species Reactivity: Human

Applications Tested: Flow cytometry (FCM)

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.1%Proclin 300

Storage

Store at 4 °C in dark.

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM

