

Anti-CD62L Mouse Monoclonal Antibody



Catalog #: UY0004

Aliases

L-selectin, LECAM-1, LAM-1, Leu-8, TQ-1

Background

Gene Name: SELL

NCBI Gene Entry: [6402](#)

UniProt Entry: [P14151](#)

Application Information

Molecular Weight: Predicted, 42 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7265

Species Reactivity: Human

Applications Tested: Flow cytometry (FCM)

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.1%Proclin 300.

Storage

Store at 4 °C in dark.

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM