

Anti-CD197 Mouse Monoclonal Antibody



Catalog #: UY0005

Aliases

BLR2; EBI1; CCR-7; CCR7; CDw197; CMKBR7; CC-CKR-7

Background

Gene Name: CCR7

NCBI Gene Entry: [1236](#)

UniProt Entry: [P32248](#)

Application Information

Molecular Weight: Predicted, 43 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7270

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

A synthesized peptide derived from human CD197

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS with 0.05% sodium azide

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1000

Flow Cytometry (FCM): 1:50-1:100

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Flow Cytometry analysis of K562 cells using CD197 antibody (green) and negative control (red).