

WB-Validated PSMB8 Lentiviral shRNA Knockdown Kit



Catalog #: V63657

Aliases

PSMB8; Proteasome 20S Subunit Beta 8; RING10; PSMB5i; D6S216E; LMP7; Multicatalytic Endopeptidase Complex Subunit C13; Really Interesting New Gene 10 Protein; Proteasome Subunit Beta Type-8; Low Molecular Mass Protein 7; Proteasome Subunit Beta 8; Proteasome Component C13; Macropain Subunit C13; EC 3.4.25.1; Beta5i; Proteasome (Prosome, Macropain) Subunit, Beta Type, 8 (Large Multifunctional Peptidase 7); Proteasome (Prosome, Macropain) Subunit, Beta Type, 8 (Large Multifunctional Protease 7); Proteasome (Prosome, Macropain) Subunit, Beta Type, 8; Large Multifunctional Peptidase 7; Proteasome Catalytic Subunit 3i; Low Molecular Weight Protein 7; Proteasome Subunit Beta 5i; Proteasome Subunit Beta-5i; Proteasome-Related Gene 7; Proteasome Subunit B5i; Protease Component C13; Proteasome Subunit Y2; D6S216; PRAAS1; ALDD; NKJO; JMP; Y2

Background

Gene Name: PSMB8

NCBI Gene Entry: [5696](#)

Storage

Store at -80 °C for one year.

Kit Components

1. WB-Validated PSMB8 shRNA lentiviral particles (4 mL)
2. Non-Target shRNA lentiviral particles (4 mL)
3. Verification Tool: KD-Validated Anti-PSMB8 Mouse mAb #63657 (5 µL)

Tested Cell Line

HT-1080

Validation Methods

RT-qPCR; Western Blotting (WB)

Shipping

Shipped with dry ice. Immediately store the product in a standard freezer at -80°C upon receipt.

Instructions For Use

The following protocol uses HT-1080 cell as an example assuming your cell culture medium is DMEM-based.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


1. Release 0.5 million HT-1080 cells into a 35 mm tissue culture dish in 2 mL of the growth medium (DMEM containing 10% FBS and 1% pen/strep). Cell density should reach 50-60% confluence the following day.
2. 24 h after cell release, pre-warm the shRNA lentiviral medium to 37°C.
3. Discard 1 mL of the original growth medium of the 35 mm dish.
4. Using a serological pipette, gently mix the lentiviral solution 3 times.
5. Carefully add 1 mL of the lentiviral solution to the well. Tip: To prevent splashing, add the solution to the dish along the wall.
6. Add a polybrene stock solution to the culture medium at a final concentration of 5 µg/mL. Gently swirl the dish to mix.
7. 48 h after cell release, without discarding the original medium, add another 1 mL of lentiviral medium directly into the dish.
8. Add an additional polybrene stock solution into the dish to obtain a final concentration of 5 µg/mL. Tip: Now, the medium in the dish should be a total of 3 mL.
9. 72 h after cell release, cells may reach confluence. Trypsinize the cells off the 35 mm dish and culture those cells in a 60 mm dish.
10. Add puromycin to the dish at a final concentration of 4 µg/mL. Tip: To assess the efficacy of puromycin selection, culture a dish of wild-type HT-1080 cells as a negative control.
11. Allow puromycin selection for 48 h. Almost all wild-type HT-1080 cells should die, while the dish infected with lentiviruses should have some remaining cells.
12. Replace the medium with regular growth medium without puromycin and allow the cells to grow to confluence before harvesting or staining.

Note: 1. This product is for research use only.

2. This product is only supplied to end users.


3. Do not use this product for commercial.

Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

RT-qPCR analysis. HT-1080 cells were infected with PSMB8-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers. $\Delta Ct (Ct_{KD} - Ct_{WT})$ was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula: $(1 - 1/2^{\Delta Ct}) \times 100\%$.



Western blotting analysis. PSMB8 protein expression in wild-type (WT) and shRNA knockdown (KD) HT-1080 cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against PSMB8 and Hsp90 α , respectively, followed by incubating with HRP-conjugated goat anti-mouse secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.